 CRESPINE® GEL+

PRILOCAINE INCLUDED

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CRESPINE® GEL+

**With 0.3 %
prilocaine hydrochloride**

Injectable implant for the treatment of osteoarthritis

**Hyaluronic Acid Implant of High Purity –
Single Injection**

Instructions for use

Description

CRESPINE®GEL+ is an absorbable intra-articular implant of high purity. It is a medical device, manufactured from hyaluronic acid of non-animal origin.

CRESPINE®GEL+ is a sterile, apyrogenic, visco-elastic, water-insoluble, biologically compatible (non-immunising, non-inflammatory, non-toxic) gel implant, manufactured from cross-linked hyaluronic acid that is obtained by fermentation. Hyaluronic acid is a naturally occurring polysaccharide in the dermal matrix of the human skin and a component of synovial fluid. Hyaluronic acid is chemically, physically and biologically identical in the tissues of all higher organisms.

Composition

1 ml CRESPINE®GEL+ contains:

Hyaluronic acid	1.0 mg
Hyaluronic acid cross-linked	14.0 mg
Sodium chloride	6.9 mg
Prilocaine hydrochloride	3.0 mg
Water for injection ad	1.0 ml

Mode of action

CRESPINE®GEL+ is injected into joints as a visco-supplement for the reduction of pain symptoms and improved functioning of the joints.

Indication and application

CRESPINE®GEL+ is indicated for pain and joint function disorders as a result of degenerative deficits as a single injection treatment course. After the injection hyaluronic acid forms a layer on the surface of the joint cartilage and improves the pain symptoms.

Treatment areas

Osteoarthritis of knee or hip.

Contraindications

CRESPINE®GEL+ must not be used in the event of:

- autoimmune illnesses
- incompatibility with gram-positive bacteria

- active inflammation or infection processes
- anticoagulant therapies
- known allergy to hyaluronic acid
- in case of known hypersensitivity against local anaesthetics of the amide type
- patients with congenital or idiopathic methemoglobinaemia or severe anaemia in children
- patients with multiple allergies should be excluded from treatment

Precautions for use

CRESPINE®GEL+ is only intended for intra-articular injection and must not be injected into blood vessels. This could lead to an occlusion of the vessels and embolism.

No clinical data is available on the administration of the product during pregnancy or lactation or on its administration to adolescents under 18 years of age.

Precautions

As with all injection treatments of the joints, the general, valid contraindications for intra-articular injection are also to be observed in treatment with **CRESPINE®GEL+**.

The skin must be disinfected prior to treatment. **As with all implants, CRESPI®GEL+ must not be used in patients with existing infections and inflammation processes in the vicinity of the implantation site.** Patients should not take any anticoagulant medication (such as aspirin), steroids or high doses of vitamin E before treatment, as these substances can lead to bleeding and a tendency to inflammation at the injection site. The syringes and needles used are considered contaminated and must therefore be destroyed in accordance with the recognised rules of medical practice.

Patients taking drugs associated with drug-induced methemoglobinemia such as sulfonamides, acetaminophen, acetanilide, aniline dyes, benzocaine, chloroquine, dapsone, naphthalene, nitrates and nitrites, nitrofurantoin, nitroglycerin, nitroprusside, pamaquine, para-aminosalicylic acid, phenacetin, phenobarbital, phenytoin, primaquine, and quinine are also at greater risk for developing methemoglobinemia; this is also true for patients with glucose-6-dehydrogenase-deficiency.

Patients with impaired renal or liver function or cardiac impairment might reach toxic concentrations of prilocaine hydrochloride earlier and should therefore be treated with caution.

CRESPINE®GEL+ should also be used with caution in patients showing symptoms of cardiac conduction disorders.

Effect on the ability to drive and use machines:

It is up to the physician to decide whether the patient is able to drive and to use machines, depending on the patient's general state of health after the intervention.

Adverse effects

Caused by the injection:

As with any other injection, patients may suffer from the following symptoms:

- temporary erythema
- slight swelling
- pain
- sensation of heat

Typically these reactions spontaneously disappear within 2 to 5 days following the injection.

Caused by the medical device:

Hypersensitivity to hyaluronic acid after the injection is reported in less than 1 %. This hypersensitivity manifests itself by extensive erythema and swelling. These reactions can occur immediately after the injection or up to 2 to 4 weeks later. The clinical data show that these hypersensitivity reactions are mild or moderate and last for a maximum of two weeks.

Caused by prilocaine hydrochloride

Possible side effects are mainly dose-dependent and inter-individually different. The following side effects have been reported from the use of prilocaine hydrochloride in higher dosages and for other indications than pain relief.

Persistent paraesthesia may be possible.

CNS manifestations like light-headedness, nervousness, apprehension, euphoria, confusion, dizziness, drowsiness, tinnitus, blurred or double vision, vomiting, sensations of heat, cold or numbness, twitching tremors, convulsions and unconsciousness are known for the use of prilocaine hydrochloride as local anaesthetic in higher concentrations.

Cardiovascular system:

Signs and symptoms of depressed cardiovascular function may commonly result from a vasovagal reaction.

Allergic reactions are characterised by cutaneous lesions, urticaria, oedema or anaphylactoid reactions. Allergic reactions as a result of hypersensitivity are rare.

Interactions

Prilocaine hydrochloride may contribute to the formation of methemoglobinaemia in patients treated with other drugs known to cause this condition (see “precautions”). The concomitant use of beta blockers may prolong the anaesthetic effect.

Application instructions

A full anamnesis must be performed before treatment, in order to rule out any contraindications. The practitioner has to advise the patient of all precautions and possible adverse effects prior to treatment. The area to be treated must be carefully prepared with an antiseptic before treatment. The syringe is removed from the blister pack, the cap covering the tip of the syringe is removed from the syringe and a needle is applied to the syringe. The needles used for implantation have to be selected according to their suitability for intra-articular injection.

Practitioners who employ local anaesthetic agents should be well versed in diagnosis and management of emergencies that may arise from their use. Resuscitative equipment, oxygen and other resuscitative drugs should be available for immediate use. To minimize the likelihood of intravascular injection, aspiration should be performed before the application. If blood is aspirated, the needle must be repositioned until no reflow of blood can be elicited by aspiration. Please note, however, that the absence of blood in the syringe does not ensure the avoidance of intravascular injection.

Advice

The injection technique is of decisive importance for the final outcome. **CRESPINE®GEL+** must only be used by authorised medical personnel in agreement with the regional legal regulations. It is possible to treat several joints at the same time. The common contraindications for intra-articular injection have to be observed.

Important information

The graduation on the syringe is an orientation aid for the user and refers to the final volume. It does not possess a measuring function, but instead only states the quantities used, relating to the nominal volume of 2 ml.

Nature and contents of container

CRESPINE®GEL+ is supplied in a 2 ml syringe with integrated Luer Lock adapter in a blister pack for single use only. The contents of the syringe are sterile (steam sterilisation). Instructions for use as well as labels with the batch number and the expiry date together with the blister pack are packed in a cardboard box. Do not use the product after its expiry date. A label is given to the patient as well in order to guarantee the traceability of the product.

Important note

The product must not be re-sterilised. It must not be mixed with other injection materials or decanted. The product must not be used if the packaging and the blister pack are damaged.

Precautions for storage

Do not freeze or expose to extreme heat.
Store at room temperature (2 °C – 25 °C).

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BioPolymer GmbH & Co. KG